


Botron High Quality Ergo-One heel grounders are made from the best products on the market. Our heel grounders have the largest cup design with better floor contact and better testing when soiled. With its dual layer construction non-marking interior and conductive outer surface is designed for long wear. The 18" conductive ribbon (snap-off tab) provides a safe static discharge to ground when tucked into your shoe. We sew our own products so we can make it any way you need. The **B7500** has stretch Velcro fastening system for comfort with a built in 1-megohm resistor, we also offer 2meg and no resistor products. **Disclaimer.** All statements of technical information are believed to be true and are based upon tests we believe to be reliable. The proper use and application for this product must be the responsibility of the user. The statements herein shall have no force or effect.



Stretch Velcro for comfort
Quick on-off attachment and removal

Installation:

1. Un-strap hook & loop strap and slip the rubber cup onto heel.
2. Pull the hook & loop tight over arch area.
3. Trim conductive ribbon to desired length and place inside of sock or shoe. (skin contact is best)

OTHER PRODUCTS:

- B7528 D-Ring heel grounder**
- B7571 Non-Marking heel grounder**
- B7525 Snap-Loc heel grounder**
- B7591 Red Hi-Visibility heel grounder**
- B7522 Toe grounder**
- Many other heel & toe grounders**

1/5" VERSES 1" CUP DESIGN

- 50% Better floor contact.**
- 50% Better chance to pass test when soiled.**
- 50% Better fit, will not slip off.**

PROPERTIES

SPECIFICATIONS

Sole Interior:	Non-Marking
Sole Exterior:	<10s
Material:	Two layer rubber
Ribbon:	18" conductive
Cup Design:	1.5" x 7.5" (others are smaller)
Charge Decay:	<0.01Sec
Thickness:	0.08-0.13mm
RTG:	with 1meg <2x10 ⁶ -1x10 ⁷
Meets:	ANSI/ESD S.20.20 & Mil 263
3/4" Stretch Hook & Loop:	Fastening system
Packing:	1 per.bag / bags of 10
Custom Sizes:	We sew our own products

RELATED TEST EQUIPMENT:

- B82250 Dual Tester with stand**
- B82450 Combo tester with stand**
- B8525 Deluxe Dual tester & stand**
- B85502 Calibration unit for testers**
- B8506 Wall test station**

Lean Stream, LLC

2153 O'Toole Avenue, Suite G
San Jose, CA, 95131
Phone: **408-894-9400**
Fax: 408-954-9989
Email: info@lean-stream.com / www.lean-stream.com

ERGO-ONE HEEL GROUNDER